



Zinc trace cuvette test, 0.02-0.8 mg/L Zn

Product #:

LCS360

Hazardous

AED Price:

Contact Hach



Hach cuvette tests: Minimum amounts of reagents, maximum safety.

Ready for use, ensures maximum safety, easy to handle, photometric evaluation system and environmentally friendly

Specifications

Digestion required:	yes, for total Zinc
Measuring range:	0.02 - 0.8 mg/L Zn
Method:	PAR
Number of tests:	24 (test includes zero solution)
Parameter:	Zinc, trace
Platform:	LCK
Special remarks:	For Total Zinc digestion use Crack Set LCW902
Storage conditions:	15 - 28 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)