



## Hydrogen peroxide pipette test 1 - 10 g/L H<sub>2</sub>O<sub>2</sub>

Product #:

LCW058

Hazardous

AED Price:

Contact Hach



For determination of Hydrogen peroxide by the Peroxomolybdate method.

---

### Specifications

Measuring range:	1 - 10 g/L H <sub>2</sub> O <sub>2</sub>
Method:	Peroxomolybdate
Number of tests:	40
Parameter:	Hydrogen peroxide
Special remarks:	additional LCW906 round cuvettes needed
Storage conditions:	2 °C - 8 °C (keep refrigerated)

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)