



## Tubing, C-flex, 6.4 mm (i.d.)

<b>Product #:</b>	<b>2327337</b>
AED Price:	Contact Hach
Available	

White, thermoplastic elastomer with excellent chemical resistance and temperature range from -80 to 230°C. Can be sterilized by gamma radiation or autoclave. Supplied in 25-foot sections.

---

### Specifications

Diameter:	5/16 in
Housing material:	C-flex
I. D.:	1/4 in
Length:	25 ft
Wall Thickness:	1/16 in