



Be Right™



Tubing, nalgene 8000, vacuum, 9.5 mm (I.D.)

Product #:	2074146
AED Price:	Contact Hach
Call for ship date	

Similar to Nalgene 8000 regular tubing, but with heavier walls. Imprinted in red for instant identification. Saves time and avoids errors which could invalidate experiments or injure laboratory personnel.

Specifications

Diameter:	7/8 in
Formula (Text):	R-3603
Housing material:	Tygon
I. D.:	3/8 in
Length:	50 ft
Wall Thickness:	1/4 in